

SpectraAnalyzer

ACCESSORY KIT 2.0 FOOD

+



ZEUTEC

EN

Accessory Kit 2.0 food

The Accessory Kit contains useful tools and a variety of sample Cups that make sample presentation and handling an easy task. The kit contains both Open and Closed Cups that help present granular and powdery samples to the SpectraAnalyzer. To optimize the filling of the Closed Cups a most helpful filling tray and a spatula are included. A brush as part of the kit simplifies the repeated cleaning of the Closes Cups. The Open Cups included are supplied for the presentation of pasteous products like chocolate or sausage meat. All accessories included in the Accessory Kit are available seperately as well.

Usage

The photostory illustrates the easy preparation of powder samples for SpectraAnalyzer analysis with the tools and sample cups included in the Accessory Kit. The fast Closed Cup filling and sample compaction procedure is illustrated. The Open Cups included in the Accessory Kit can also be placed into the sample tray and filled with any pasteous products.



Specifications

Scope of Supply

Transport case	Brush (10306)
Two Closed Cups (10206)	Spatula with wooden handle (10300)
One Open Cup (10201)	Sample loader assy (10303)

Order information

Accessory Kit 2.0 food	100-B9001-1
------------------------	-------------

ZEUTEC Opto-Elektronik GmbH

Friedrich-Voß-Straße 11, 24768 Rendsburg

T (+49) 4331 - 136650

E moreinfo@zeutec.de

www.spectraalyzer.com

ZEUTEC