


Sample presentation accessories



SpectraAlyzer standard drawer models:



Closed Cup

for the analysis of powders: e.g. milk powder, starch, ground coffee and flour.



Open Cup

for the analysis of pastes: e.g. butter, cream cheese and minced meat.



Open Cup deep

for the analysis of olive paste, olive cake and minced meat.



Slide Cup

for the analysis of products with low water content: e.g. oil and glycerine.



Italian Cup

for the analysis of chocolate mass and for clinical analysis of human faeces.



British Cup

for the analysis of liquid eggs.



Viscous Cup

for the analysis of viscous samples: e.g. honey, ketchup and condensed milk.



Cheese Cup

for the analysis of hard cheese

SpectraAlyzer rotating drawer models:



Rotating Cup

for the analysis of granular samples:
e.g. grains, corn or pellets.



Rotating Cup deep

for the analysis of granular samples
and fluffy materials: e.g. cotton or wool.



Rotating Open Cup flat meat

for the analysis of sticky materials:
e.g. minced meat.



Rotating Open Cup deep meat

for the analysis of sticky materials:
e.g. minced meat.

SpectraAlyzer Top Window models:



Cup Top Window

for the analysis of powders: e.g. milk
powder, starch, ground coffee, flour
and olive paste.



Liquid Cell Assy Top Window

for the analysis of liquid samples: e.g.
oil, glycerine, marinades and ketchup.

Accessories list

| Product name | Item no. |
|---------------------------------|-------------|
| 1. Closed Cup | 10206 |
| 2. Open Cup | 10201 |
| 3. Open Cup deep | 10202 |
| 4. Slide Cup | 10203 |
| 5. Italian Cup | 10209 |
| 6. British Cup | 10217 |
| 7. Viscous Cup | 10210 |
| 8. Cheese Cup | 10260 |
| 9. Rotating Cup | 10242 |
| 10. Rotating Cup deep | 10241 |
| 11. Rotating Open Cup flat meat | 10252 |
| 12. Rotating Open Cup deep meat | 10251 |
| 13. Cup Top Window | 400-BH014-1 |
| 14. Liquid Cell Assy Top Window | 410-BH001-1 |

ZEUTEC Opto-Elektronik GmbH

Friedrich-Voß-Straße 11
24768 Rendsburg
Germany

(+49) 4331 - 136650
moreinfo@zeutec.de
www.spectraalyzer.com

ZEUTEC