

## SpectraAlyzer Printer

### Description



**SpectraAlyzer - Printer** is a portable thermal printer that offers an innovative alternative to impact based systems. Very light weight and compact size, easy automatic paper loading, end of paper and black mark sensors.

**SpectraAlyzer - Printer** is powered using an external self-detecting power supply, with RS232 serial interface.

### Technical data

<b>Printing method</b>	Thermal line printing
<b>Resolution</b>	203 x 203 dpi
<b>Printing speed</b>	> 50 mm/sec
<b>Columns</b>	24
<b>Character set</b>	ASCII standard, International
<b>Characters/line</b>	32 (12x24); 42 (9x24)
<b>Paper width</b>	58 mm
<b>Roll dimension</b>	57.5 mm ± 1mm
<b>Roll thickness</b>	60 mm
<b>Interface</b>	RS232
<b>Power</b>	100 – 240 VAC
<b>Reliability MTBF</b>	100 km
<b>Safety</b>	EN60950+A1+A2+A3+A4
<b>Dimension</b>	146 x 88 x 65 mm
<b>Weight</b>	370 g (without paper roll)

ZEUTEC Opto-Elektronik GmbH

Friedrich-Voß-Straße 11  
24768 Rendsburg  
Germany

(+49) 4331 - 136650  
moreinfo@zeutec.de  
www.spectraalyzer.com

**ZEUTEC**