

# Quality check and control

Pet food / wet



**ZEUTEC**

 made  
in  
Germany



# One **SpectraAlyzer MEAT** for all parameters!

## FAST, EASY, ACCURATE & MODULAR

Analytical parameters and performance are:

Parameter	Range [%]	Error (+ / -) *	Repeatability
Moisure	23.0 - 30.0	0.25	< 0.002
Moisure	75.0 - 86.0	0.50	< 0.002
Protein, low	1.85 - 10.79	0.35	0.01
Protein, all	1.8 - 37.3	0.64	0.01
Fat, all	1.7 - 21.5	0.51	0.01
Fat, low	1.7 - 8.3	0.39	0.01
Ash	1.0 - 10.2	0.43	< 0.002

\* **RMSEP (Root Mean Squared Error of Prediction)**

**ZEUTEC Opto-Elektronik GmbH**

Friedrich-Voß-Straße 11  
24768 Rendsburg  
Germany

(+49) 4331 - 136650  
moreinfo@zeutec.de  
www.spectraalyzer.com

**ZEUTEC**