

Quality check and control

Scalded / parboiled sausage and ham



ZEUTEC

 made
in
Germany



One **SpectraAlyzer MEAT** for all parameters!

FAST, EASY, ACCURATE & MODULAR

Analytical parameters and performance are:

Parameter	Range [%]	Error (+ / -) *	Repeatability
Moisure	53.0 - 70.0	0.47	< 0.05
Protein	10.4 - 16.8	0.33	< 0.05
Fat	13.7 - 30.0	0.56	< 0.05
BEFFE **	8.8 - 14.0	0.43	< 0.05
Ash	2.0 - 4.1	0.08	< 0.005

* **RMSEP (Root Mean Squared Error of Prediction)**

** **BEFFE % = Protein % - (foreign Protein % + connective tissue [Collagen] Protein %)**

Connective tissue Protein % = Hydroxyprolin content x 8

ZEUTEC Opto-Elektronik GmbH

Friedrich-Voß-Straße 11
24768 Rendsburg
Germany

(+49) 4331 - 136650
moreinfo@zeutec.de
www.spectraalyzer.com

ZEUTEC