

Quality check and control

Cured sausage and ham



ZEUTEC

 **made
in
Germany**



One **SpectraAnalyzer MEAT** for all parameters!

FAST, EASY, ACCURATE & MODULAR

Analytical parameters and performance are:

Parameter	Range [%]	Error (+ / -) *	Repeatability
Moisture	22.2 - 58.1	0.77	< 0.05
Protein	13.7 - 31.6	0.58	< 0.05
Fat	11.9 - 48.6	0.62	< 0.05
BEFFE **	11.9 - 28.5	0.58	< 0.01
Ash	5.0 - 6.0	0.06	< 0.01

* **RMSEP (Root Mean Squared Error of Prediction)**

** **BEFFE % = Protein % - (foreign Protein % + connective tissue [Collagen] Protein %)**

Connective tissue Protein % = Hydroxyprolin content x 8

ZEUTEC Opto-Elektronik GmbH

Friedrich-Voß-Straße 11
24768 Rendsburg
Germany

(+49) 4331 - 136650
moreinfo@zeutec.de
www.spectraalyzer.com

ZEUTEC