

ZEUTEC SpectraAlyzer GRAIN NEO Analysenparameter

Überblick Getreide, Leguminosen, Ölsaaten

Produkt	Wasser (%)	Protein (%)	Hektoliter (kg/L)	Öl (%)	Stärke (%)	Nassgluten (%)	Sedimentation (ml)	Alveo W	Kornhärte (HI)	Asche (%)	Faserstoffe (%)
Gerste	•	•	•		•						
Weizen	•	•	•		•	•	•		•		
Weizenmehl	•	•			•	•	•	•		•	
Hartweizen/Durum	•	•	•			•					
Hirse	•	•	•								
Mais	•	•	•	•	•						
Hafer	•	•	•								
Reis	•	•									
Reisschalen	•	•		•	•						•
Sojabohnen	•	•	•	•							•
Kichererbsen	•	•									•
Linsen	•	•									
Raps	•	•		•							
Senfsaat	•			•							
Roggen	•	•	•								

