

ZEUTEC SpectraAlyzer GRAIN NEO Analytical Parameters

Overview for grains, pulses and oilseeds

Product	Moisture (%)	Protein (%)	Hectoliter (kg/L)	Oil (%)	Starch (%)	Wet Gluten (%)	Sedimentation (ml)	Alveo W	Hardness (HI)	Ash (%)	Fiber (%)
Barley	•	•	•		•						
Wheat	•	•	•		•	•	•		•		
Wheat Flour	•	•			•	•	•	•		•	
Durum	•	•	•			•					
Sorghum	•	•	•								
Corn	•	•	•	•	•						
Oat	•	•	•								
Rice	•	•									
Rice Bran	•	•		•	•						•
Soybean	•	•	•	•							•
Chickpea	•	•									•
Lentil	•	•									
Canola	•	•		•							
Mustard Seed	•			•							
Rye	•	•	•								