

Quality and improved Yield



ZEUTEC

 made
in
Germany



One **SpectraAlyzer DAIRY** for all parameters!

FAST, EASY, ACCURATE & MODULAR

Designed as a modular system, the SpectraAlyzer DAIRY solution presents the analytical results of these major quality parameters within 45 seconds.

Application examples

- **Milk** - Fat, protein, lactose, total solids, casein, titrable acidity, SNF
- **Milk powder** - Moisture, fat, protein, lactose, ash, acidity, SNF
- **Pre-mix and ice cream** - Fat, lactose, SNF, MSNF, protein
- **Butter and Spreads** - Moisture, protein, fat, salt
- **Yoghurt** - Fat, total solids, acidity and protein
- **Cheese** - Fat, moisture, casein, solids not fat
- **Khoya** - Moisture, fat, protein, total solids, ash, total energy
- **Cream** - Moisture, fat, protein, total solids, ash, total energy
- **Whey** - Fat, solids, protein, lactose
- **Dried whey** - Moisture, fat, titrable acidity, total solids, ash
- **Paneer** - Moisture, fat, protein, total solids, salt

ZEUTEC Opto-Elektronik (India) Private Limited

Sabi Cottage, 681, Kavery Layout,
HA Farm Post Dasarahalli,
Bangalore – 560024, Karnataka,
India

T: (+91) 7003752371
E: marketing@zeutec.com

ZEUTEC GmbH

Friedrich-Voß-Straße 11
24768 Rendsburg
Germany

T: (+49) 4331 - 136650
E: moreinfo@zeutec.de